ILLINOIS REGION 3 PROTOCOLS

COMBITUBE AIRWAY (ALS)

I. Indications:

- A. Cardiac/respiratory arrest
- B. Respiratory arrest comatose with no gag reflex

II. <u>Contraindications:</u>

- A. Patient with an intact gag reflex
- B. Patient under age 16 and/or under 5 feet tall
- C. Patient with known esophageal disease
- D. Patient with a history of esophageal trauma/or ingestion of caustic substance
- E. Patient with a tracheostomy or laryngectomy
- F. Patient with a foreign body obstruction in the trachea

III. Potential Complications:

A. Damage to the proximal cuff from broken teeth or dentures

IV. Precautions:

A. DO NOT remove the Combitube in the field unless the patient's gag reflex returns or the patient has been endotracheally intubated. Remove only upon the order of Medical Control.

V. <u>Equipment:</u>

- A. Combitube airway
- B. 100 ml syringe
- C. 15 ml syringe
- D. Water based lubricant
- E. Suction equipment

VI. Procedure:

- A. Hyperventilate with 100% oxygen
- B. Assemble equipment, check cuffs for leaks
- C. Place head in neutral position
- D. Insert device in midline using jaw lift

PROCEDURE CONTINUED

- E. Insert until black rings are at the teeth
- F. Inflate the pharyngeal cuff with 100 ml air (blue pilot balloon)
- G. Inflate the distal cuff with 15 ml air (white pilot balloon)
- H. Ventilate through the blue tube
- I. Auscultate lung sounds and over the epigastrium
- J. If bilateral breath sounds are auscultated,, continue ventilation through blue tube
- K. If there is absent chest rise, no lung sounds and gurgling heard over the epigastrium, begin ventilating through the clear tube labeled #2.
- L. Confirm placement. Auscultate lung sounds and over the epigastrium. Observe for chest rise.

Blessing Hospital, EMS Medical Director	Passavant Hospital, EMS Medical Director
	Tussavani Hospital, Elvis Medical Bricetor
St Johns Hospital, EMS Medical Director	Mamarial Madical EMC Madical Director
	Memorial Medical, EMS Medical Director